

Focus on your working style.



ZEISS ZEN core

Your Microscope Software for User-specific Repetitive Workflows in Biomedical Labs

zeiss.com/zen-core



Seeing beyond

Your Microscope Software for User-specific Repetitive Workflows in Biomedical Labs

› In Brief

› The Advantages

› Packages and Features

› The Applications



› Technology and Details

› Service

You decide what your microscope software looks like.

With ZEN core you can adapt the user interface exactly to your application.

- Configure your ZEN core only with necessary functions.
- Define and combine workbenches easily to create workflows, so called jobs.
- Jobs guide your operators through a defined flow of consecutive tasks ensuring data repeatability and ease of use.
- Combine Job Mode with user management functions to assign specific imaging tasks to certain users and user groups.
- Opt for the GxP or Multi Channel toolkit, when traceability or fluorescence imaging of biological samples is important to you.



The ZEN core user interface provides both a bright and a dark mode to meet the needs of different users and their preferred working environment.

Simpler. More Intelligent. More Integrated.

- › In Brief
- › **The Advantages**
- › Packages and Features
- › The Applications
- › Technology and Details
- › Service

Easy to Configure. Easy to Use.



ZEN core gives you the benefit of an adaptive user interface tailored to your needs of a biomedical lab. The easy-to-follow GUI configuration accommodates tasks of all kinds and any complexity. Define visible functions as well as dedicated working procedures. ZEN core also offers you configurable user management, so you can specify users and user roles. Whatever their level of experience, operators will learn the software quickly.

Define workflows. No workarounds.

With ZEN core you can set up individual workflows – from image acquisition, processing, analysis to final reports. Use the built-in automated image acquisition routines, such as HDR or Best Image to benefit from the consistency of a repeatable workflow. With application-specific modules you enhance your microscope to get typical workflows of your biomedical lab simply done.

You need GxP compliance? Go for it.

Are you working in pharmaceutical or medicine production or in another regulated environment? Do you require FDA 21 CFR Part 11 conformity? Then, count on full traceability along your microscopy workflow with the GxP-toolkit. It provides save and conform means for reliable quality results.



Expand Your Possibilities

- › In Brief
- › **The Advantages**
- › Packages and Features
- › The Applications
- › Technology and Details
- › Service

User Management Designed to Assure Data Repeatability and Integrity




Supervisor
Experienced User



Operator
Routine User



Administrator
Lab Manager



Sample configuration of different user roles


Expand Your Possibilities

- › In Brief
- › **The Advantages**
- › Packages and Features
- › The Applications
- › Technology and Details
- › Service

Keep doing your workflows

Don't be afraid of repetitive tasks. Define a Job and setup your imaging workflow. Define, lock, or hide microscope parameters according to your needs. Only the weekend can stop your operators to happily get reproducible results. Day after day.

Workflow Example



Your Simple and Affordable ZEN Packages for Laboratory Routine

- › In Brief
- › The Advantages
- › **Packages and Features**
- › The Applications
- › Technology and Details
- › Service

Instruments	Manual Light Microscope	Motorized Light Microscope	SEM/FIB-SEM
Base Software	ZEN starter	ZEN core	
Acquisition Packages	Base	Motorized	Smart Acquisition
Toolkit Packages	2D Toolkit	3D Toolkit	Vision4D Toolkit *
	Connect Toolkit	AI Toolkit	GxP Toolkit
System / Driver / Special	Coded Mic	Linkam	


* Sold as a stand-alone software

ZEISS ZEN core at Work

- › In Brief
- › The Advantages
- › Packages and Features
- › **The Applications**
- › Technology and Details
- › Service

GxP Toolkit

The GxP toolkit meets the requirements of regulated industries, such as pharmaceutical companies that often require FDA 21 CFR Part 11 compliance. It allows you to trace complete workflows through seamlessly integrated microscopy hardware and software. No matter if you need to characterize the particle size of powders, starch, or other substances, if you check for the presence of contaminants, or if you monitor the microbial load of a production facility (gram staining), with ZEN core, every workflow can be made GxP compliant.



Time	User	Category	Description	Comment
09.04.2019 16:56:00	Supervisor	ExitMode	Login	
09.04.2019 16:56:01	Supervisor	EnterMode	Manage Templates	
09.04.2019 16:56:18	Supervisor	EnterMode	Create Job Template	
09.04.2019 16:56:24	Supervisor	ExitMode	Create Job Template	
09.04.2019 16:56:25	Supervisor	EnterMode	Manage Templates	
09.04.2019 16:56:50	Supervisor	EnterMode	Create Job Template	
09.04.2019 16:57:20	Supervisor	Save	The following item has saved: Routine 1	
09.04.2019 16:57:36	Supervisor	ExitMode	Create Job Template	
09.04.2019 16:57:37	Supervisor	EnterMode	Manage Templates	
09.04.2019 16:58:08	Supervisor	SetStatus	An electronic signature was made on job template 'Routine 1' using cert-ID '1'-, 'D1'- On Supervisor@PC134117003192468178001640816 Commitment made: 'http://url.letslog.org/01903/v1.2.24/ProofOfCreation' (created by me.garten@zeiss.com) 'GxP'	
09.04.2019 16:58:08	Supervisor	SetStatus	Status of template 'Routine 1' changed from 'Draft' to 'Actor' (Electronic Signature)	
09.04.2019 16:58:08	Supervisor	ExitMode	Manage Templates	
09.04.2019 16:58:30	Supervisor	ExitMode	Free Examination	
09.04.2019 16:58:30	Supervisor	Logout	User 'Supervisor' has logged off	
09.04.2019 16:58:32	Operator	Login	User 'Operator' has logged in	
09.04.2019 16:58:32	Operator	ExitMode	Login	
09.04.2019 16:58:34	Operator	EnterMode	Run Job Template	
09.04.2019 16:58:46	Operator	EnterMode	Run one time	
09.04.2019 16:58:46	Operator	Execute	Execute the template 'Routine 1' with status: 'Actor' In task '1' 'ZD Acquisition' the following parameter has changed: In tool Camera the following parameter has changed: ExposureTime has changed from '100' to '20'	because image was too bright
09.04.2019 16:59:07	Operator	Execute	In task '1' 'Interactive Measurements' the following parameters have changed:	

Audit Trail of all user activities

The GxP toolkit helps you to:

- Guarantee reproducible results and secured processes to prevent regulatory consequences
- Automatically document your microscopy examinations to make your audit safe
- Provide backup / disaster recovery functionality, and archive / database capability

Use a variety of tools and functionalities in combination with the required qualification and validation activities to maintain GxP compliance for your images, tables, and reports:

- User management
- Traceable audit trail
- Release procedure of workflows
- Electronic signatures, incl. countersign functionality
- Checksum protection of process-critical data

ZEISS ZEN core at Work


- › In Brief
- › The Advantages
- › Packages and Features
- › **The Applications**
- › Technology and Details
- › Service

ZEN Data Storage and Data Explorer

Do you acquire many microscope images or data? Do you need to collaborate or consult with internal or external partners? Remember how long it took, when you have searched for a certain data set that was generated 2 years ago? ZEN Data Storage with Data Explorer is your solution for smart storage and data access. You create a central image and data repository that is accessible from anywhere and anytime.

Work efficiently and modern while you:

- Use an intelligent, central data management solution
- Effortlessly share images, data, reports and workflows
- Access data with any mobile device or PC, even on the other end of the globe
- Easily browse, view and annotate images
- Use role based permissions when working with partners



ZEISS ZEN core at Work

- › In Brief
- › The Advantages
- › Packages and Features
- › **The Applications**
- › Technology and Details
- › Service

ZEN Multi Channel Acquisition*


Fluorescence microscopy is a standard technique in biomedical disciplines and helps understanding structure and dynamics within cells, tissue, and other biological samples. Being able to automatically acquire and overlay fluorescence signals in individual channels is key for optimal visualization and data analysis. Prominent examples for fluorescence imaging include the characterization of the spatial distribution of proteins, structural changes, cell-cell interactions, and colocalization experiments.

ZEN Multi Channel Acquisition allows for:


- Easy & intuitive definition of your experiment using Smart Setup
- The combination of fluorescence and transmitted light channels
- Channel specific camera and light source settings

ZEN Multi Channel Acquisition can also be combined with the GxP toolkit if you work in regulated environments or need to ensure data authenticity. In Job mode, complete workflows including image acquisition, analysis, and reports can be realized.

* Part of Base Acquisition Toolkit



2D Multi-Channel Acquisition workbench allows to acquire images with multiple (fluorescent) channels in one go with automated switching of illumination and camera mode settings per channel.



Smart Setup allows to easily define the channels of the acquisition.

Your Flexible Choice of Components

- › In Brief
- › The Advantages
- › Packages and Features
- › The Applications
- › **Technology and Details**
- › Service

General Base Functionality			
General Functions	ZEN starter	ZEN core	Description
General Workbenches	●	●	General image acquisition and image analysis functionality via workbenches
Job Mode		●	Configurable mode to guide operators through a pre-defined workflow
Basic Measurement	●	●	Basic functions for interactive measurement
Measurement	●	●	Extended functionality of interactive measurement
Online Measurement		●	Measurement of areas and non-area values in the live image
Topo View		●	3D visualization of topography images and hightmap representation
Reporting	●	●	Creation of reports containing detailed information
Report Template Creator (MS Word® Add-In)	●	●	Software to create report templates using MS Word®
ZEN Connect Entry	●	●	Interactive acquisition and contextual display of images in single- and multi-instrument workflows
CAD Import		●	Import CAD data to create overlay images
ZEN Data Storage Client	●	●	Seamlessly connect to ZEN Data Storage Server for central storage of documents and templates
LM-specific			
	ZEN starter	ZEN core	Description
Microscope Control		●	Control of light microscopes and devices via software
Camera Control	●	●	Control of cameras from software
Panorama	●	●	Manual acquisition of high resolution images; automated acquisition and stitching functionality on coded and non-coded stages
Manual Extended Focus	●	●	Manually acquire multiple images at different focus positions and combine them to an image with a greater depth of field
EM-specific			
	ZEN starter	ZEN core	Description
Microscope Control		●	Control of scanning electron microscopes and devices via software

● included ○ optional

Your Flexible Choice of Components

- › In Brief
- › The Advantages
- › Packages and Features
- › The Applications
- › **Technology and Details**
- › Service

Acquisition Toolkits		ZEN starter	ZEN core	Description
Base Acquisition	Multichannel Acquisition			<ul style="list-style-type: none"> ■ Enable routine applications: ■ Simultaneous acquisition of multiple (fluorescence) channels ■ Acquisition of images at defined intervals
	Time Series	○	○	
Motorized Acquisition	Motorized Extended Focus			Enable motorized applications on your microscope: <ul style="list-style-type: none"> ■ Automatically acquire multiple images at different focus positions and combine them to create an image with a greater depth of field ■ Determine the focus position of your sample ■ Record exact and high-resolved images of large samples by automatically scanning pre-defined areas
	Autofocus		○	
	Tiles & Positions			
Smart Acquisition	Guided Acquisition		○	Automatically identify interesting positions in overview images and create detail scans of respective positions
Linkam	Linkam (Hardware Control)		○	Interface to control Linkam heating and cooling stages
Coded Microscope	Coded microscope	○		Read out of coded microscope components

Toolkit		ZEN starter	ZEN core	Description
GxP	GxP		○	Ensure traceability and accountability of workflows (precondition for 21 CFR Part 11 compliance)
2D Toolkit	Image Analysis		○	Toolkit for 2D image analysis by creation of automatic measurement programs. Advanced processing functions included
	Advanced Processing			
3D Toolkit	3Dxl		○	Visualization of 3D data in ZEN core
	Advanced Processing			
AI Toolkit	ZEN Intellesis Segmentation			Complete package for AI application including integrated training interfaces: <ul style="list-style-type: none"> ■ Automated image segmentation based on machine learning algorithms ■ Automated object classification of segmented and analyzed images based on machine learning algorithms ■ AI-based denoising of images using noise-2-void algorithms
	ZEN Intellesis Object Classification		○	
	ZEN Intellesis AI denoising			
Connect Toolkit	ZEN Connect			Complete Connect package to extend the ZEN Connect functionality. Includes L-marker calibration for correlative LM to EM workflows and functionality for 3D data. Import 3 rd -party microscopy images and metadata into ZEN core using Bio-Formats included.
	ZEN Connect 2D Add-on			
	ZEN Connect 3D Add-on		○	
	Third-party Import			

● included ○ optional

Your Flexible Choice of Components

- › In Brief
- › The Advantages
- › Packages and Features
- › The Applications
- › **Technology and Details**
- › Service

System Requirements	
Minimum System Requirements	
Processor	Intel® CPU with 6 cores, 3.3 GHz clock frequency, e.g. Intel® Core® i5-10600 CPU
Memory	16 GB RAM
Graphics adapter	Support for resolutions of at least 1920 × 1200 pixels, 32 Bit true color, DirectX® 12.0 (or higher), e.g. integrated Intel® UHD Graphics 630 or dedicated graphic card NVIDIA Quadro P1000 4GB
Hard drive	1× 512 GB SSD NVMe and 1× 8 TB HDD SATA 6Gb/s
USB interface	2× USB 2.0 and 2× USB 3.0
Monitor	24" TFT display, vertical resolution 1200 pixels at least
Operating system	Microsoft® Windows® 10 x64 Professional or Enterprise, no special customer adapted versions.

Recommended System Requirements	
Processor	Intel® Xeon® CPU with 4 cores, 3.8 GHz clock frequency, e.g. Intel® Xeon® Gold 5222
Memory	64 GB RAM
Graphics adapter	Support for resolutions of 3840 × 2160 pixels ("4K"), 32 Bit true color, DirectX® 12.0 (or higher), OpenGL® 4.5 (or higher), e.g. dedicated graphics card NVIDIA Quadro RTX4000 8 GB.
Hard drive	1× 512 GB SSD NVMe and 2× 6 TB HDD SATA 6Gb/s (configured as 6 TB RAID 1)
USB interface	4× USB 2.0 and 2× USB 3.0
PCIe interface	2 available PCI Express x8 Gen2 slots for camera interface cards
Monitor	27" TFT display, vertical resolution 1440 pixels at least
Operating system	Microsoft® Windows® 10 x64 Professional or Enterprise, no special customer adapted versions.

ZEISS Service – Your Partner at All Times

Your microscope system from ZEISS is one of your most important tools. For over 170 years, the ZEISS brand and our experience have stood for reliable equipment with a long life in the field of microscopy. You can count on superior service and support - before and after installation. Our skilled ZEISS service team makes sure that your microscope is always ready for use.

- › In Brief
- › The Advantages
- › Packages and Features
- › The Applications
- › Technology and Details
- › **Service**

Procurement

- Lab Planning & Construction Site Management
- Site Inspection & Environmental Analysis
- GMP-Qualification IQ/OQ
- Installation & Handover
- IT Integration Support
- Startup Training

Operation

- Predictive Service Remote Monitoring
- Inspection & Preventive Maintenance
 - Software Maintenance Agreements
 - Operation & Application Training
 - Expert Phone & Remote Support
 - Protect Service Agreements
 - Metrological Calibration
 - Instrument Relocation
 - Consumables
 - Repairs

New Investment

- Decommissioning
- Trade In

Retrofit

- Customized Engineering
- Upgrades & Modernization
- Customized Workflows via APEER



Please note: Availability of services depends on product line and location

>> www.zeiss.com/microservice



Carl Zeiss Microscopy GmbH

07745 Jena, Germany
microscopy@zeiss.com
www.zeiss.com/zen-core